

INVISIBLE IN CANAL (IIC)

Our new IIC - Invisible-in-canal - is a clever, all-in-one design


INVISIBLE IN CANAL



Despite its small size, the IIC comes with outstanding Widex technology.

Not almost invisible, but invisible

The IIC model is placed deep into the ear canal so no part of it can be seen.



A CIC model is placed in the external part of the ear canal. Our IIC model however, is placed beyond the second bend of the ear canal.

Because it is tailor-made for each individual hearing aid user with our unique CAMISHATM system, the IIC fits perfectly.

Outstanding features of the IIC include:

- **Perfect fit:** each IIC hearing aid is custom made to ensure a perfect fit.
- **Better localisation:** tight fit means users can better determine where sounds are coming from
- **All-purpose use:** can be worn in all listening situations
- **Durability:** all components are contained in one unit making it extremely durable
- **More comfort:** no joined parts mean there are no hard edges
- **Clever design:** Attractive and stylish, all the components are contained in one unit.

And of course, as a Widex product, it benefits from the Widex Sound - natural sound that doesn't just help users hear speech or loud sounds, but all the sounds around them.



The IIC (Invisible In Canal) model is now available, call us for more details on [01202 861522](tel:01202861522)

