

Accessories, Materials,
Supplies & Tools



MICR((()))**SONIC**[®]

CLEARLY THE SOUND CHOICE



Accessories, Materials, Supplies & Tools

Welcome to Microsonic's Accessories, Materials, Supplies & Tools Catalog, created to help make your job as a hearing health care professional a little easier.



Microsonic facility in Ambridge, PA

For Your Convenience...

Eight well-organized sections are packed with hundreds of unique and essential items...many new...from around the world. The product photos and detailed descriptions are designed to help you find and order the product you need quickly and easily.

Global Supplier

Commitment to highest quality, on-time delivery, and competitive pricing have made Microsonic a global supplier of hearing health industry from its US and European locations. We export our products to over 60 countries in four continents. Microsonic has become one of the largest hearing aid tubing manufacturer in the world.

Quality

With ISO 9001 and I3485 certifications, our facilities are committed to deliver highest quality hearing health products. The great majority of the products in this catalog are manufactured in US and Europe.

Innovative

Microsonic has developed many products and processes over the decades. Our contemporary R&D facility houses numerous state-of-the-art equipment and experienced personnel to bring you novel products.

Your ongoing support of Microsonic is genuinely appreciated. We have more than four decades of experience as a full-service earmold lab, and are one of the industry's leading suppliers of top-quality impression material and hundreds of other hearing health care-related products.

When you have a technical question, or aren't sure which product will best suit your client's needs...ask us. Chances are our knowledgeable staff will be able to make a recommendation on the spot.

Please let us know if you need additional copies – we'll be pleased to send them. And please let us know of new items that you find that are beneficial to your practice, so we can evaluate them for inclusion in our Accessories and Supplies Catalog.

Five Easy Ways to Order

Your order can be placed any of five convenient ways:

Register at online store
to access
special pricing for professionals

1 Call us toll-free at 1-877-376-7139

2 Fax your order to us at 1-724-266-9482

3 Mail your order to us at
2960 Duss Ave. Ambridge, PA 15003

4 E-Mail accessorydepartment@microsonic-inc.com

5 Shop online at <http://store.microsonic-inc.com>

Use the Five-Point Check List

To insure quick processing of your order, and to avoid errors, enter the following information for each item on your order before you transmit it to us:

✓ Item Number(s) from catalog	✓ Name of item(s)	✓ Quantity of each item	✓ Size(s) where applicable	✓ Color(s) where applicable
--	----------------------------	----------------------------------	-------------------------------------	--------------------------------------

Sizes and colors appear in the Online Store as separate product numbers – be sure to select the correct item!

Pricing

A current Accessory Price List may be obtained by contacting us. Prices are subject to change without prior notice. Additional updated product and/or pricing information may be sent to you periodically.

Invoicing Terms

Net 30 days upon credit approval. Accessory items are not discounted. An invoice will accompany your shipment and a statement will be mailed on or about the 28th of each month. Past due accounts may be subject to credit hold and late fee.

Shipping and Handling

We ship most orders via UPS Ground, though small orders under \$10.00 may be sent First Class Mail. Please call for actual shipping charges. In general, ground shipping and handling will be charged as follows:

Order Amount	Shipping/Handling
\$0-\$50.00	\$7.75
\$50.01-\$100.00	\$8.75
\$100.01-\$250.00	\$9.75
\$250.01-\$500.00	\$11.95
\$500.01-\$1000.00	\$12.95
\$1000.01+	Please call

These rates are for UPS Ground in the 48 Continental US only. COD shipments are typically charged normal weight-based shipping rates plus the UPS COD charge. Shipping rates are subject to change without prior notice.

When Speed is Needed


Express deliveries can be arranged to most locations. All shipping costs are billed for this service. Same-day UPS shipment is available ONLY for orders received by 2:00PM Eastern Time. Same-day FedEx or shipments must be placed by 12 Noon Eastern Time.

Returns

Please promptly inspect the product upon delivery. We will replace or credit your account for any **NON-CONFORMING PRODUCT** within warranty guidelines. ALL ACCESSORY ITEMS REQUIRE A RETURN AUTHORIZATION NUMBER THAT CAN BE OBTAINED BY CALLING 1-877-376-7139. Returns may be subject to a 5% restocking fee.

Growing for you Online

New products are continually added to bring you the best choices possible. Be sure to check our Online Store (store.microsonic-inc.com) regularly for new items that may not appear in our print catalogs. Should you not find the product you are looking for, please let us know. We are here to serve you, and may be able to add it to our selection!



SafeTouch, templet, MegaPress, Megapress+, Mega-Sil, Quikfloaters, EZ-Gel, and Advanced Design Horn (A.D.H.) are trademarks of Microsonic, Incorporated. Tube-Lock Plus and Advanced Design Horn (A.D.H.) are patented by Microsonic, Incorporated.

Microsonic Incorporated does not accept liability arising out of the use, the misuse or the inability to use any accessory, tool or product described herein, or any accessory, tool or product subsequently offered for sale by Microsonic, Incorporated. The user assumes all risk and liability whatsoever in connection with these accessories, tools and products.

Table of Contents

Impression Supplies, Materials & Tools

Micro-Sil™ Silicone Impression Material 8

Quikfloaters™ Silicone Compound.....8

Mega-Sil™ Silicone Impression Material 8

Empty Solo Container 8

Mega-Sil™ Solos 8

templet™ Impression Material 9

instaPRO™ Instant Hearing Protection Material 9

Megapress™ Impression Material 9

Silasoft™ 9

Stainless Steel Syringe 9

Microsonic Universal Syringe – Standard 10

Audalin Acrylic Impression Material 10

Kot-N-Blox 10

Oto-Blocks 10

CIC Oto-Blocks 10

Blunt End Tweezers 10

Bite Blocks™ 11

Diagnostic Otoscope 11

Advanced Slimline Earlight 11

Cerumen Management Tips 11

Replacement Parts for Advanced Slimline Earlight 11

Oto-Probe Magnifying Adapter 11

Bionix® Lighted Ear Curette™ with Magnification 12

Microsonic Universal Syringe CIC 12

Purple Knight Gloves 12

Heine Mini Fiber Optic Otoscope 12

Earmold Tubing, Supplies & Accessories

Coiled Tube 14

N.A.E.L. Standard Pre-Bent Tubing 14

Tube-Lock Assembly 14

Tube-Lock Plus (U.S. Patent 5,488,205) 15

Invisi-Lock 15

Moisture Free Tubes 15

Libby Horn Tubes 16

Continuous Flow Adapter (C.F.A.™) Tubes 16

Silicone Stethoscope Tubing 16

Tubing Cement 16

Tubing Puller 16

Loctite® Super Bonder 495 16

Tubing Expander 17

Connectors-Elbows 17

Receiver Attachment Connectors 17

Pilot Communicator 17

Tubing Scissors 17

Synth-a-Flex II Polishing Kit 17

Acoustic Modification & Venting Accessories

Adcobuild – For Major Build-up 20

Hard Adco Addon – For Minor Build-up 20

Soft Adco Addon Build – For Minor Build-up 20

Knowles Acoustic Dampers 20

Select-a-Vent (S.A.V.) Trees 21

ITE Personal Listener 21

Testing Accessories

Stethoscope 24

Non-Toxic Test Tips 24

E-A-RLink® Foam Tips for Insert Earphones 24

E.A.R® Demonstration Tips 24

Comply™ Canal Tips 24

Comply™ Snap Tips 25

Comply™ Contoured Earpieces 25

Comply™ Snap Ear Piece 25

Ear Level Testing Stethoscope 25

Safesound Earphone Covers 25

Tools, Cleaning & Finishing Products

Tubing Air Blower 28

Jewelers' Screwdrivers (Set of 3) 28

X-Acto® Knife 28

Spring-loaded Needle-Nose Pliers 28

Improved Pin Vise Reamer Tool 28

3" Polishing Wheel 28

3" Buffing Wheel 28

Polishing Bar Compound 29

Grinding Stones 29

Oto-Slik™ 29

Wax-Away Tool Kit 29

Wax Removal Tools 29

Vent Cleaning Tool 29

Jiffy Bath Contact Cleaner 30

Wondertool® ITE Battery Extractor ... 30

Disposable Vent Brush 30

Oto Floss 30

Safetouch™ Wipes 30

Ultrasonic Solution 30

Alligator Forceps 30

5" Stainless Steel Hemostat 31

Waxman 2000 31

Contact Renu Cleaner 31

Ear Care, Hearing Aid & Earmold Maintenance

Comply™ Soft Wraps 34

Ad-Hear™ Cerumen Guards 34

Ad-Hear™ Plus Cerumen Guards 34

SafeTouch™ Disposable Currettes 34

SafeTouch™ Earmold Cleaning Kit 34

SafeTouch™ Bulb Syringes 34

SafeTouch™ Drying Kit for Hearing Aids 35

SafeTouch™ Perfectdry® Electric Drying System 35

SafeTouch™ Cleaning Spray 35

SafeTouch™ Cleansing Wipes 35

SafeTouch™ EZ-Gel™ 35

SafeTouch™ Starter Set 35

SafeTouch™ Parents' School Kit 36

Super Dri-Aid™ Kit – Box of 6 36

Audiologist's Choice Anti-Itch Cream ... 36

Dri At Nite 36

Miracell 36

Otoferm 36

Tech-Care Ear Wax Removal Drops ... 36

PerfectRite RIC Cleaning System

Main Function 38

Commercial Adaptation 40

SafeTouch™ PerfectRite® RIC Cleaning System Configurations41

SafeTouch™ PerfectRite® Adapter for Universal RIC/RITE Eardomes41

SafeTouch™ PerfectRite® Cleaning Cartridges41

SafeTouch™ PerfectRite® Custom Adapter Kit41

SafeTouch™ PerfectRite® Desiccants41

Miscellaneous & Educational

Swim Plug Lanyards 43

Mini Battery Tester 43

Activair Button Cell Tester 43

Digital Hearing Aid Battery Tester/Holder 43

ETY-PLUG Hi-Fi Generic Attenuators 43

Ear Planes™ 43

Communication Collar Clip 43

Telephone Pads™ 44

Hearing Aid Sweat Bands 44

Super Seals™ 44

Ear Band-It™ 44

Perfect Earmold Impressions DVD 44

Mounted Demonstration Ear 45

Resealable Plastic Bags 45

Instruction Ears 45

Premium Mailing Boxes 45

See-Thru Plastic Pouches 45

Leatherette Pouch 45

Microsonic Mailing Box 45

Earmold Order Book 46

E-A-R™ Classic Earplugs 46

E-A-R™ Express Pod Plugs 46

Gator Clips 46

Gator Clips II 46

Quick Reference Guide to Earmold Selection 46

Replacement Handles 47

Universal Hearing Screener™47

ER Attenuators 47

Self-Inking Stamp 47

Microsonic Custom Earmold Manual – 9th Edition 47

Study Guide To Accompany Microsonic's Custom Earmold Manual 47



What's new in Microsonic Accessory Catalog?

Completely renewed SafeTouch™ hygiene product line for variety of in-ear devices

Our laboratories work with hearing professionals to develop unique cleaning solutions to help preserve your hearing systems. We tested each product at every stage to ensure our customers the most reliable and effective results.



Novel impression material formulas for Hearing Professionals

Multi-purpose **MegaPress+**™ and high precision fast set **templet**™ impression materials are available in two-part 50ml cartridges:



also:

- ✓ New Clamshell Packaging for tubing
- ✓ The Earmold Manual 9th Edition
- ✓ PerfectRite® RIC Cleaning System
- ✓ Fluxx Audio instaPRO™ instant hearing protection material
- ✓ New Quikfloaters™ Cartridge Floatable Swim Plug Material

Impression Supplies,
Materials & Tools





MICRO-SIL™ SILICONE IMPRESSION MATERIAL

A premium silicone impression material, Micro-Sil™ is very easy to handle and highly accurate. The catalyst blends smoothly with the base material in 20 to 30 seconds and cures in 5 to 8 minutes. Micro-Sil™ impressions will never sag or collapse and cannot be damaged by handling or shipping.

Micro-Sil™ is available in a bulk kit yielding approximately 58 impressions. Available in purple. Fully cured Shore-A: 25±3. Imported from Germany.

- #100 Micro-Sil™ Bulk Kit with Tube Catalyst
- #101 Micro-Sil™ Base only
- #102 Tube Catalyst only

QUIKFLOATERS™ SILICONE COMPOUND

Quikfloaters™ is a superior two-part vinyl polysiloxane impression compound specifically formulated for custom in-office swim molds. Packaged in bulk or pre-measured, pocket-sized individual units. Quikfloaters™ is available in red or blue. Please specify color when ordering. Tear-resistant, floatable. Cures in 5 minutes. Fully cured bulk kit Shore-A: 22±3; Cartridge system Shore-A: 35, curing time 10 minutes

- #104 Quikfloaters™ Bulk Kit (makes approx. 10 pairs)
- #105 Quikfloaters™ Solo Individual Pack (specify red or blue)
- #106 Quikfloaters™ Sealant
- #112 Quikfloaters™ Impression Gun
- #173 Quikfloaters™ Cartridge System (Red & Blue – 4 cartridges/ea.)

MEGA-SIL™ SILICONE IMPRESSION MATERIAL

Mega-Sil™ is our best two-part vinyl polysiloxane system, specifically compounded for custom earmold impressions for both ITE instruments and CIC fittings. An excellent soft material, Mega-Sil™ cures in approximately 5 minutes and produces non-greasy impressions which are extremely accurate and dimensionally stable. Mega-Sil™ is available in yellow, peach or blue. Please specify color when ordering. Fully cured Shore-A: 30±3. Made in USA.

- #107 1.3 lbs. 35 impressions (specify yellow, peach or blue)
- #108 2 lbs. 54 impressions (specify yellow, peach or blue)
- #109 4 lbs. 108 impressions (specify yellow, peach or blue)

EMPTY SOLO CONTAINER

For packaging your own impression material.

- #110 Empty Solo Container

MEGA-SIL™ SOLOS

Pocket-sized Mega-Sil™ Solos are individual pre-measured units of Mega-Sil™ premium vinyl polysiloxane. No measuring – just open and blend for a perfect mix every time. Fast-curing, highly accurate and shrink-resistant Mega-Sil™ Solos are excellent for standard-depth and CIC impressions. Comes in a handy reusable pack. Available in yellow, peach or blue. Shore-A: 30±3

- Mega-Sil™ Solos Individual Pack (specify yellow, peach or blue)
- #111B Yellow
- #111C Blue
- #111E Peach



templet™ IMPRESSION MATERIAL

The new templet™ is a polydimethylsiloxane impression material designed specifically to take multiple impressions in a short period of time, and achieve excellent details. No measuring or mixing needed. The material automatically mixes as it moves through the tip of the syringe. One cartridge will make 4-5 impressions. Cures in 2 minutes. Shore-A: 38. Available in distinct turquoise color.

- #160 templet™ Impression Gun
- #165 templet™ Impression Cartridge
- #162 templet™ Nozzle (specify small, medium or large)

instaPRO™ INSTANT HEARING PROTECTION MATERIAL

Very comfortable floatable silicone material for direct fitting of individual hearing protection from FLUX Audio. instaPRO™ is designed to provide extreme tensile strength and durability. Cures in approx. 5 minutes. Final hardness Shore A 35. Available in Red, blue, orange, green, yellow, white, and black.

- #112 instaPRO™ Impression Gun
- #174 instaPRO™ Cartridge material
- #114 instaPRO™ Nozzle (specify medium or large)
- #174K instaPRO™ Kit (2x red, 2x blue, 1x orange, 1x green, 1x yellow, 32 medium nozzles)

MEGAPRESS+™ IMPRESSION MATERIAL

MegaPress+™ is a premium multi-purpose impression material, perfect for taking multiple impressions in a short period of time. Available in disposable cartridges. No measuring or mixing needed. One cartridge will make 4-5 impressions. Cures in approximately 3 minutes. Highly accurate and shrink-resistant MegaPress™ impressions are excellent for standard depth and CIC impressions. Available in distinct orange color. Fully cured Shore-A: 39±3

- #160 MegaPress™ Impression Gun
- #161 MegaPress™ Impression Cartridge
- #162 MegaPress™ Nozzle (specify small, medium or large)

SILASOFT™


Silasoft™ is a vinyl polysiloxane impression material formula with remarkable softness. It offers the same smooth blending qualities, easy insertion and excellent memory as other vinyl addition silicone systems, but features a lower Shore-A rating of 28±3 (fully cured). This additional softness makes it ideal for sensitive ears and CIC fittings. Cures in 5 minutes. Each kit yields 40 impressions. Made in USA.

- #115 Silasoft™ Impression Kit
- #115A Silasoft™ Solos (box of 24)

STAINLESS STEEL SYRINGE

This is the ultimate in professional earmold impression syringes. Meticulously crafted of medical grade stainless steel, it will last a lifetime. Perfect to use with both powder/liquid and silicone materials, it is self-cleaning of even uncured residue. Ships with matching nontip base. Specify standard or C.I.C. nozzle.

- #121 Stainless Steel Syringe
- #122 Standard Nozzle
- #123 C.I.C. Nozzle



128R

MICROSONIC UNIVERSAL SYRINGE – STANDARD

A patented Microsonic exclusive. Our comfortable, high quality plastic syringe with one-piece molded gripper helps reduce hand fatigue. Use with all impression materials.

- #128 Microsonic Universal Syringe (Pink)
- #128R Replacement “O” Ring



116

AUDALIN ACRYLIC IMPRESSION MATERIAL

This traditional impression material is available in a bulk package containing material for about 45 impressions, or in a kit of 24 pre-measured units which automatically provide the perfect proportion for a good impression. Cures in 10 minutes. Cured Shore-A: 26±3.

- #116 Audalin Bulk Kit
- #117 Audalin Powder only
- #118 Audalin Liquid only
- #120 Empty containers to package your own Audalin (specify powder or liquid type and quantity)




132

KOT-N-BLOXS

Made of cotton, Kot-N-Blox are pre-threaded and pre-tied. Available in 4 convenient sizes. 50 per package. Specify size.

- #132 Mini
- #133 Small
- #134 Medium
- #135 Large




136

OTO-BLOCKS

Super-soft Oto-Blocks foam canal blocks are pre-threaded for easy, safe use and removal. Available in 4 universal sizes. 50 per package. Specify size.

- #136 Small 7mm x 8mm high
- #137 Medium 9mm x 8mm high
- #138 Large 11mm x 8mm high
- #139 X-Large 13mm x 8mm high



140

CIC OTO-BLOCKS

The shorter CIC Oto-Blocks are specifically designed to help protect the delicate inner canal wall and eardrum during the CIC impression procedure. 3 sizes. 50 per package. Specify size.

- #140 Small 7mm x 6mm high
- #141 Medium 9mm x 6mm high
- #142 Large 11mm x 6mm high



143

BLUNT END TWEEZERS

Rounded safety ends for inserting cotton or foam canal blocks prior to impression taking.

- #143 Blunt End Tweezers



144

BITE BLOCKS™

For open-jaw impressing. Package of 100

- #144 Bite Blocks



145

DIAGNOSTIC OTOSCOPE

Imported from Germany, our Heine Mini-Set medical grade otoscope is extremely lightweight and a pleasure to use. Uses two AA batteries. Additional disposable specula can be ordered separately in packages of 50.

- #145 Diagnostic Oscope Mini-Set
- #146 50 Replacement Disposable Specula (specify 2.5 or 4mm)
- #147 Replacement Light Bulb



146-2



149

ADVANCED SLIMLINE EARLIGHT

This top quality earlight features a high-intensity krypton bulb for excellent visibility. Uses AAA size batteries, included.

- #149 Advanced Slimline Earlight




150

CERUMEN MANAGEMENT TIPS

Specially designed to aid in the management of cerumen, these tips fit our Advanced Slimline Earlight.

- #150 Cerumen Management Tips



151

REPLACEMENT PARTS FOR ADVANCED SLIMLINE EARLIGHT

- #151 Replacement Tip for Advanced Slimline Earlight
- #152 Replacement Bulb for Advanced Slimline Earlight



153

OTO-PROBE MAGNIFYING ADAPTER

Just snap on to the Advanced Slimline Earlight for a magnified view of the ear canal and eardrum.

- #153 OTO-Probe Magnifying Adapter

154




BIONIX® LIGHTED EAR CURETTE™ WITH MAGNIFICATION

This distinctive system features single-use curettes that snap easily into a brilliant light source for superb viewing in the canal. A magnifier is provided as well. Available in 5 tip styles. 50 per box, including 1 light source and 1 magnifier. Light Source maybe purchased separately.

- | | |
|----------------------------|---|
| #154 Lighted Angleloop® | #154D Lighted InfantScoop® |
| #154A Lighted Microloop® | #154E Lighted Variety Pack – All 5 Styles |
| #154B Lighted WaveCurette® | #154F Light Source |
| #154C Lighted FlexLoop® | |

163



MICROSONIC UNIVERSAL SYRINGE CIC

A patented Microsonic exclusive. Our comfortable, high quality plastic syringe with one-piece molded gripper helps reduce hand fatigue. Use with all impression materials.

- #163 Microsonic Syringe for CIC Impressions (Blue)

164



PURPLE KNIGHT GLOVES

Nitrile gloves protect hands from damage, and may be safely used to mix silicone impression materials. 100 per box. Purple only. Five sizes available. Please specify size (XS, S, M, L, or XL).

- #164 Purple Knight Gloves

167



HEINE MINI FIBER OPTIC OTOSCOPE

Same convenient size as the popular miniature diagnostic otoscope, but with higher-visibility fiber optic head. Uses standard 146-2 and 146-4 specula.


- #167 Heine Mini Fiber Optic Otoscope
#168 Replacement Bulb for Heine Mini Fiber Optic Otoscope



168

Earmold Tubing,
Supplies & Accessories

MICR(SONIC®



200

COILED TUBE

This compact tubing stretches and recoils as needed. For use with communication earmolds. Extends up to 36".

#200 Coiled Tube (specify right or left)



201

N.A.E.L. STANDARD PRE-BENT TUBING

Highest quality earmold tubing available in N.A.E.L. (National Association of Earmold Laboratories) Standards. Please specify color when ordering: clear, pink or brown.

#201 Single Bend Quilled Tubes. Package of 50 and 100 (specify size)
#202 Double-Bend Quilled Tubes (specify quantity, right or left & size)
#203 Bulk Tubing. Available by the foot. (specify length & size)

Standard N.A.E.L. Tubing (In inches)


Size	Type	Outside Dia.	Inside Dia.
12	Standard	.125	.085
13	Medium	.122	.076
13	Thick	.130	.076
13	Double Wall	.142	.076



202



203



204

TUBE-LOCK ASSEMBLY

An exclusive Microsonic design, the Tube-Lock coupling system is an excellent way to attach tubing to soft earmolds. When inserted into a soft earmold, a tiny gold-coated flanged ring on the quilled tubing fastens into the material forming a secure seal, eliminating the need for cementing. Ideal for vinyl, silicone and soft silicone earmolds. For easy installation, use the Tube-Lock Applicator Tool. Package of 10, 50, and 100 Tube-Lock Assemblies. Available in clear, pink or brown.

#204 #13 Thick Tube-Lock Assembly
#204A Gold Lock for #13 Thick
#204D #13 Thick Tubes with Silver Lock
#204E Silver Lock #13 Thick and Moisture Free
#205 #13 Double Wall Tube-Lock Assembly
#205A Brass Lock for #13 Double Wall
#206 Tube-Lock Applicator Tool



204D



204E



205A

206



207

207B




208


TUBE-LOCK PLUS (U.S. Patent 5,488,205)

The Tube-Lock Plus is a coupling system with patented, pressure-fit retaining flanges, designed to work with soft earmolds. The ridges hold the tube firmly in place without the need for cements or adhesives, yet allow tubing replacement as needed. Package of 10, 50, and 100 Tube-Lock Assemblies. Available in clear, pink or brown.

#207 Tube-Lock Plus
#207B Tube-Lock Plus – Lock Only
#208 Tube-Lock Plus Double Wall



207E




207F


INVISI-LOCK

Microsonic's Invisi-Lock is similar to the Tube-Lock Plus but without the surface flange that limits insertion depth. It enhances cosmetic appeal while retaining its gripping capabilities. Package of 10, 50, and 100 Tube-Lock Assemblies. Available in clear, pink or brown.


#207E Invisilock – Lock Only
#207F Invisilock Tubing Assembly



209



210



211



212



213



214, 215



216



217

MOISTURE FREE TUBES

Condensation of moisture in earmolds and tubing is a problem for some clients. Moisture Free is a denser tubing material that greatly reduces and often eliminates this moisture build-up. **N.A.E.L. (National Association of Earmold Laboratories) Standards sizes.** Package of 10, 50, and 100 Tube-Lock Assemblies. Available in clear, pink or brown.

#209 #13 Thick Single Bend (specify quantity)
#210 #13 Thick Double Bend (specify quantity & right or left)
#211 #13 Thick Tube-Lock Assembly (specify quantity)
#212 #13 Thick Tygon (bulk) (by the foot, specify length)
#213 3mm or 4mm Libby Horn (specify quantity & size)
#214 #13 Medium C.F.A.™ (specify quantity)
#215 #13 Thick C.F.A.™ (specify quantity)
#216 #13 Medium Single Bend (specify quantity)
#217 Moisture Free Tube-Lock Plus (specify quantity)




221-223

LIBBY HORN TUBES

The Libby Horn is a smooth, tapered, one-piece sound tube of internal stepped bore design based on Killion acoustic horn principles. Available in clear, tint or brown. *Specify quantity.*

- #221 3mm Regular Libby Horn
- #222 3mm Thick Libby Horn
- #223 4mm Libby Horn



CONTINUOUS FLOW ADAPTER (C.F.A.™) TUBES

A replacement C.F.A.™ tube and connector provide the convenience of a snap-in connection. Tubing can never crimp or change in diameter. Available in clear, tint or brown. *Specify quantity.*

- #224 #13 Medium C.F.A.™ Tubes
- #225 #13 Thick C.F.A.™ Tubes
- #226 #13 Double Wall C.F.A.™ Tubes
- #227 Reverse Curve C.F.A.™ Tubes
- #247 C.F.A.™ Seating Rings



228

SILICONE STETHOSCOPE TUBING

A rugged, non-shrinking tubing for use with stethoscope adapter earmolds or anesthesiologist's earpiece adapters. Available by the foot.

- #228 Silicone Stethoscope Tubing (*specify length*)



230

231

TUBING CEMENT

For strong bond when attaching tubing to an earmold. **CAUTION: ALWAYS use with adequate ventilation.**

- #230 Tubing Cement for Acrylic & Semi-Soft Molds
- #231 Tubing Cement for Synth-A-Flex II Molds



232

TUBING PULLER

For easy insertion of tubing through any earmold, use the handy nylon loop Tubing Puller.

- #232 Tubing Puller



233

LOCTITE® SUPER BONDER 495

Cyanoacrylate adhesive for those tough gluing jobs. 10 fluid oz. (3ml). (*Does not glue Medi-Sil II, M2000 or Polyethylene materials.*)

- #233 Loctite® Super Bonder 495



234

TUBING EXPANDER

This tool temporarily expands the tubing to make it easier to slide onto any type of accessory, attachment or connector.

- #234 Tubing Expander



235

236

CONNECTORS-ELBOWS

- #235 Long Shank/Polyethylene Clear
- #236 Short Shank/Polyethylene Clear



243

245

RECEIVER ATTACHMENT CONNECTORS

- #243 Male/Polyethylene Clear
- #245 Female/Polyethylene Clear



248A

248

PILOT COMMUNICATOR

Holds a pilot's microphone in place. Typically used with our #28 earmold.

- #248 Crimped Bar
- #248A Straight Bar



250

TUBING SCISSORS

Heavy Duty stainless steel scissors allow firm control in cutting earmold tubing to the desired length.

- #250 Tubing Scissors



542

SYNTH-A-FLEX II POLISHING KIT

Quickly smooths and polishes Synth-A-Flex earmolds after adjustments are made. Includes a polishing cloth. 2 fluid ounces. **CAUTION: Use with adequate ventilation and rubber gloves. Wash earmold thoroughly after polishing.**

- #542 Synth-A-Flex II Polishing Kit

Acoustic Modification & Venting Accessories

MICR(SONIC®



ADCOBUILD – FOR MAJOR BUILD-UP

Simple to use and cures in minutes. This powder and liquid combination is mixed, applied directly to the hearing aid or earmold, and cured in the ear for a fit no other build-up product can provide. Contains enough material for many in-office repairs. Note: ADCOSHEEN should be applied over cured ADCOBUILD to obtain a seal and prevent discoloration.

- #300 Adocbuild (1/2oz. Liquid, 1oz. Powder)
- #301 Adcosheen (1/20z. Liquid)
- #302 Adcosheen Thinner (1/2oz. Liquid)



HARD ADCO ADDON – FOR MINOR BUILD-UP

Brush on a build-up for acrylic molds or ITE cases. Use when a hard finish is desired. Changes the size, not the shape of the earmold or hearing aid. Dry completely before inserting in the ear.

- #303 Hard Adco Addon Build and Thinner



SOFT ADCO ADDON BUILD – FOR MINOR BUILD-UP

Brush on in coats when a soft, flexible finish is desired. Will change the size, not the shape of the earmold or hearing aid. Dry completely before inserting in the ear.


- #304 Soft Adco Addon Build – For Minor Build-Up



KNOWLES ACOUSTIC DAMPERS

Used in earmold tubing for damping. Available in six color-coded values of acoustic resistance. Supplied in packages of 10 (one ohm value per package). Specify ohms.

- #318 680 ohms (white) Acoustic Dampers
- #319 1000 ohms (brown) Acoustic Dampers
- #320 1500 ohms (green) Acoustic Dampers
- #321 2200 ohms (red) Acoustic Dampers
- #322 3300 ohms (orange) Acoustic Dampers
- #323 4700 ohms (yellow) Acoustic Dampers



SELECT-A-VENT (S.A.V.) TREES

- #331 S.A.V. Trees (specify clear or tint)
- #332 Mini S.A.V. Tree

Each S.A.V. Tree contains one of each vent size:
#1(.156") #4(.062")
#2(.125") #5(.031")
#3(.095") and one solid plug #6




ITE PERSONAL LISTENER

A convenient listening device which fits over the canal tip of the ITE aids. Connects to your personal earmold.

- #414 ITE Personal Listener


Testing Accessories



STETHOSCOPE

Durable steel stethoscope used to test all types of hearing aids. Does not require a personal earpiece for the listener.


- #400 Stethoscope
- #401 Spare Bell Tips for Stethoscope



NON-TOXIC TEST TIPS


Non-Toxic Test Tips are available in sizes 1-6 (small - large) and fit either the right or left ear. Each test tip comes with elbow and tubing attached.

- #404 Non-Toxic Test Tips



E-A-RLINK® FOAM TIPS FOR INSERT EARPHONES

- #406A 3A Standard Pkg. of 50
- #406B 3B Small Pkg. of 50
- #406C 3C Large Pkg. of 24
- #406D Infant Eartips 3.4mm
- #406E Infant Eartips 4.0mm



E.A.R.® DEMONSTRATION TIPS

Disposable foam tips that compress to a small diameter, then expand to the size of the ear canal. Used for demonstration and testing. 24 tips with tubing attached per package.

- #407 E.A.R.® Demonstration Tips



COMPLY™ CANAL TIPS

Use Comply™ Canal Tips when a loaner is necessary for broken ITEs, for hearing aid evaluations, out of office calls and demonstrations, or trial use of a hearing aid.


- #408 Professional Fitting System
- #409 Professional Fitting System Expanded
- #412 Twelve-Packs (Specify size)
- #413 Sound Tube Adapters



COMPLY™ SNAP TIPS

Convenient for testing; perfect for ears that cannot tolerate an earmold. For use only with the Comply™ Snap Earpiece. (#416) Available in 3, 6 and 9mm depths and three widths; small, medium and large. *Please specify sizes.*

- #415 Comply™ Snap Tips



COMPLY™ CONTOURED EARPIECES

Quickly and easily attach Comply™ tips to BTE hearing aids.

- #416L Left
- #416R Right



COMPLY™ SNAP EAR PIECE

Essential adapter for use with Comply™ Snap Tips. One size.


- #416 Comply™ Snap Ear Piece



EAR LEVEL TESTING STETHOSCOPE

Instantly determine whether a hearing aid has weak or distorted performance. Invaluable for parents to check a child's hearing aid before sending to school each day.

- #420 Ear Level Testing Stethoscope



SAFESOUND EARPHONE COVERS

Economical, disposable earphone covers provide a clean surface for each patient. Box of 100.

- #426 Safesound Earphone Covers

Tools, Cleaning &
Finishing Products

MICR(SONIC®



229

TUBING AIR BLOWER

Removes moisture and foreign matter from earmold tubing.

#229 Tubing Air Blower



500

JEWELERS' SCREWDRIVERS (Set of 3)

These miniature screwdrivers feature swivel heads, aluminum handles and steel blades in the three most commonly used sizes: .055", .070" and .100". Overall length is 4".

#500 Jewelers' Screwdrivers



502

X-ACTO® KNIFE

No. 1 Knife has a 4-7/8" aluminum handle with a No. 11 blade for delicate cutting of lightweight material.

#502A No. 11 Blade 5/pkg



503

SPRING-LOADED NEEDLE-NOSE PLIERS

Needle Nose pliers with convenient spring action for ease of use.

#503 Spring-loaded Needle-Nose Pliers



504

IMPROVED PIN VISE REAMER TOOL

Easy-to-grip reamer tool comes with two collets to accommodate various drill sizes. Four drill bits are included:

#504 Improved Pin Vise Reamer Tool

#35 (.110") #49 (.070")
#43 (.089") #56 (.043")



519

3" POLISHING WHEEL

This 3" cloth wheel is used to smooth adjustments or remove scratches. *For use with acrylic and semi-soft molds only.*

#519 3" Polishing Wheel



520

3" BUFFING WHEEL

This 3" cloth wheel adds a beautiful shine to earmolds. *For use with acrylic and semi-soft molds only.*

#520 3" Buffing Wheel



521

POLISHING BAR COMPOUND

Apply to a dry polishing wheel to smooth adjustments or polish out scratches. *For use with acrylic and semi-soft molds only.*

#521 Polishing Bar Compound



523A

523

524

525

GRINDING STONES

Description, For Use With Materials:

#523 very coarse M2000, Medi-Sil II, M-2000 Plus, Neon Lites (cone shaped)
#523A very coarse M2000, Medi-Sil II, M2000 Plus, Neon Lites (bulb shaped)
#524 coarse Synth-A-Flex II, Vinylflex II, Polyethylene
#525 smooth Acrylic, Alerite



529

528

OTO-SLIK™

This grease-free, non-staining, super slippery-liquid comfortably eases insertion of earmolds and ITEs into the ear. Packaged in a handy vial that fits easily in pocket or purse.

#528 Oto-Slik™ (Individual)
#529 Oto-Slik™ (Display Case of 12)



531

WAX-AWAY TOOL KIT

Three stainless steel mini-reamers for cleaning wax and dirt from vents, earmold/receiver tubing and microphone ports.

#531 Wax-Away Tool Kit




532

WAX REMOVAL TOOLS

Designed to safely clean accumulated wax from the receiver tube and vent of an ITE aid or from the sound bore of an earmold. Set includes brush, loop and pick. *Each tool may also be purchased separately.*

#532 3-Piece Wax Removal Tool Set
#533 Wax Loop only
#534 Wax Brush only
#535 Wax Pick only
#536 Magnetic Wax Loop



534

533

536

537

535

VENT CLEANING TOOL

Economical plastic tool for cleaning wax from earmold or earshell vents.

#537 Vent Cleaning Tool

538

**JIFFY BATH CONTACT CLEANER**

This cleaner provides dependable, long-lasting protection for instrument contacts and controls. Cleans and lubricates. Sprays away noise due to dust, dirt and corrosion.

#538 Jiffy Bath Contact Cleaner

539

**WONDERTOOL® ITE BATTERY EXTRACTOR**

For stuck battery removal, even when battery door has broken. Convenient design helps extract batteries from all hearing aids.

#539 Wondertool® ITE Battery Extractor

540

**DISPOSABLE VENT BRUSH**

A quick way to clean vents through earmolds or ITE aids.

2.5mm (pkg of 10)
5.0mm (pkg of 10)

#540 Disposable Vent Brush

541

**OTO FLOSS**

For cleaning the tubing that connects the hearing aid to the earmold. Rigid end penetrates the tubing, then pulls a flexible cleaning element through. Package of 30.

#541 Oto Floss

545



546

**SAFETOUCH™ WIPES**

Disinfectant towelettes designed for use on hearing aids, earmolds, headphones, probe tips – any surface or object that needs to be cleaned or disinfected.

#545 SafeTouch Wipes™ (160)
#546 SafeTouch Wipes™ Porta-Pack (40)

548

**ULTRASONIC SOLUTION**

For use with Ultrasonic Cleaning Machine, this solution safely cleans earmolds and other small objects. Eight ounce concentrate makes 16 quarts. (1/2 oz. concentrate to one qt. of water.)

#548 Ultrasonic Solution

571

**ALLIGATOR FORCEPS**

Useful for extracting wayward cotton or foam blocks if the strings breaks.

#571 Alligator Forceps

573

**5" STAINLESS STEEL HEMOSTAT**

Great for holding tubing and and other small parts. Think of it as a third hand.

#573 5" Stainless Steel Hemostat

576

**WAXMAN 2000**

This handy vacuum device helps remove moisture from hearing aids and earmolds. Two tip sizes are included.

#576 Waxman 2000

675

**CONTACT RENU CLEANER**

Safe, efficient spray solution for cleaning the surface of batteries. Renews contact point, extends usable life of battery. Do not inhale or get in eyes.

#660 Contact Renu Cleaner

Ear Care, and Hearing Aid & Earmold Maintenance

MICR(())SONIC®

336



COMPLY™ SOFT WRAPS

Wrap around an ITE or ITC hearing aid to assure a tight, yet comfortable seal. These adhesive-backed, disposable wraps come 10 per package.

#336 Comply™ Soft Wraps

337



AD-HEAR™ CERUMEN GUARDS

A simple way to prevent wax-related hearing aid failure. When applied, a spongy material covers the speaker tube to keep wax out, yet allow amplified sound to pass transparently through to the eardrum. Easy application tab. Secure adhesive sticks to all earmold materials. 10 per pack.

#337 Ad-Hear™ (9mm)
#338 Ad-Hear™ (7mm)
#339 Ad-Hear™ (5mm)

340



AD-HEAR™ PLUS CERUMEN GUARDS

Same as Ad-Hear™ Guards with a larger spongy protective area and application tab. 10 per pack.

#340 Ad-Hear™ Plus (4mm)
#341 Ad-Hear™ Plus (6mm)
#342 Ad-Hear™ Plus (8mm)

700



SAFETOUCH™ DISPOSABLE CURETTES

Soft, strong and convenient to use, these low cost disposable currettes feature a unique, flexible basket-like tip to gather wax easily and quickly. Convenient flip-top dispenser contains 50 large and 50 small currettes.

#700 SafeTouch™ Disposable Curette
#701S SafeTouch™ Disposable Currettes Refill Small
#701L SafeTouch™ Disposable Currettes Refill Large



762A

SAFETOUCH™ EARMOLD CLEANING KIT

SafeTouch™ Earmold Cleaning Kit is ideal for disinfection. Cleaning tablets are specifically designed to dissolve earwax. Simply place your earmold and a cleaning tablet into the cup, and add water. Your earmolds are clean and free of earwax in approx 30 minutes.

#762 SafeTouch™ Earmold Cleaning Kit, Display case of 12
#762A SafeTouch™ Earmold Cleaning Kit
#763 SafeTouch™ Earmold Cleaning Tablet Refills, Box of 20

704



705



SAFETOUCH™ BULB SYRINGES

Safety collar design permits easy flushing of ear while preventing accidental insertion deep in the ear canal. Available in 1 oz. and 2 oz. sizes.

#704 SafeTouch™ Safety Syringe, 1oz.
#705 SafeTouch™ Safety Syringe, 2oz.



766A

SAFETOUCH™ DRYING KIT FOR HEARING AIDS

NEW

Perspiration and moisture cause oxidation on electronic components. Using Drying Kit on daily basis is a cost effective routine to protect your expensive Hearing Devices. Can be reactivated in microwave or oven.

#766 SafeTouch™ Drying Kit, Display case of 12
#766A SafeTouch™ Drying Kit



765

SAFETOUCH™ PERFECTDRY® ELECTRIC DRYING SYSTEM

NEW

PerfectDry® will eliminate all humidity and condensation inside your hearing aids in approx. 2 hours. Secure drying cycle with timer (automatic switch-off). Max temperature is 104° F, Max air ventilation without condensation.

#765 SafeTouch™ PerfectDry® Electric Drying System

761A



SAFETOUCH™ CLEANING SPRAY

NEW

SafeTouch™ Cleaning Spray is designed to be a safe and effective way to thoroughly clean your hearing aids. The spray works by dissolving earwax on the surface of earmold or ITE unit. To ensure cleanliness, the soft brush is packed in an individual blister pack. Bactericidal & Fungicidal, fresh mint scent. 30ml.

#761 SafeTouch™ Cleaning Spray 30ml, Display case of 12
#761A SafeTouch™ Cleaning Spray 30ml

764A



SAFETOUCH™ CLEANSING WIPES

NEW

SafeTouch™ Cleansing wipes are designed for occasional maintenance. They are individually packed for quick and easy use. Ideal for travel, they also comply with the new aviation regulations. Bactericidal & Fungicidal, fresh mint scent. Contains 30 individual packs.

#764 SafeTouch™ Cleansing Wipes, Display case of 12
#764A SafeTouch™ Cleansing Wipes, Box of 30 individually packed wipes

760A



SAFETOUCH™ EZ-GEL™

NEW

For natural soothing and calming. The specially developed SafeTouch™ EZ-Gel™ makes the fitting and removing of the Hearing Aid, Earmold or Hearing Protector easier. The EZ-Gel allows a better tolerance of the Hearing Aid and prevents irritations and skin problems. It is formulated to minimize the risks of allergies. 5ml Tube.

#760 SafeTouch™ EZ-Gel™ 5ml, Display case of 20
#760A SafeTouch™ EZ-Gel™ 5ml

768



SAFETOUCH™ STARTER SET

NEW

The SafeTouch™ Starter Set is perfect for on-the-go users. Easy and convenient storage; this re-usable bag contains everything you will need to keep your hearing devices clean & dry.

Contents: 1 Soft bag, 10 Cleansing Wipes, 1 Dry Kit, 1 Cleaning Spray, 1 EZ-Gel™ Cleaning tools, User Guide.

#768 SafeTouch™ Starter Set



921

SAFETOUCH™ PARENTS' SCHOOL KIT NEW

A convenient kit to enable parents to check the performance of their child's hearing aid before or after the school day.

Contains SafeTouch™ Starter Kit in addition to stethoscope, tubing air blower, and battery tester. Ideal for school nurses' use, also.

#921 Parents' School Kit



622-6

SUPER DRI-AID™ KIT – Box of 6

Renewable silica gel desiccant kit that removes moisture from hearing aids.
CAUTION: Remove batteries from hearing aids before using kit..

#622 Super Dri-Aid™ Kit – Individual
#622-6 Super Dri-Aid™ Kit – Box Of 6



630

AUDIOLOGIST'S CHOICE ANTI-ITCH CREAM

Hydrocortisone formula helps to relieve itchy ears with no sting, no pain. Maximum strength allowed without a prescription. For external use only.

#630 Audiologist's Choice Anti-itch Cream



665

DRI AT NITE

Disposable containers for the overnight removal of moisture from hearing aids. Consumer discards when indicator turns pink. Remove hearing aid battery before using. (Box of 12)

#665 Dri At Nite



527

MIRACELL

A concentrated blend of botanicals to help relieve pain, redness, discomfort and swelling in the ears. 0.5 fl. oz. bottle

#527 Miracell



530

OTOFERM

This cream eases insertion and helps improve comfort of earmolds, while providing an effective acoustic seal.

#530 Otoferm



702

TECH-CARE EAR WAX REMOVAL DROPS

Safe, fast-acting formula to use as an aid to soften, loosen and remove ear wax.

#702 Tech-Care Ear Wax Removal Drops Only



PerfectRite®

The first and only
RIC cleaning solution
to eliminate earwax

PerfectRite®

PerfectRite® is a patented Medical Device, with CE marking, for the cleaning and drying of RITE/RIC Hearing Aids as well as ITE.

Main function of the PerfectRite®

PerfectRite® has been designed to enable patients maintaining RIC/RITE/ITE hearing aids on their own. The user only needs to place the Hearing Aids in PerfectRite®. Within only 2 hours, the Hearing Aids will be cleaned, disinfected and dried. Your patients will benefit the full performance of their devices every day.




"Forget earwax inconvenience, dirty Domes and clogged Receivers."

1 All of our tests have been made with natural earwax Earwax on various devices



2 Start operation _____ Cleaning _____ After 2 hours _____ Satisfaction



3 Earwax removed



Commercial adaptation

PerfectRite®: The Custom Solution for RIC/RITE/ITE hearing aid user

You can customize PerfectRite® for any size of universal eardomes, as well as custom made RIC/RITE earmolds and ITEs with custom adapter fabricating kit.

The PerfectRite® Adapters are designed to fit all receivers manufactured by major brands.



PerfectRite®: Satisfaction & Loyalty of your customers

You can assure your patients for the best performance, they can easily maintain their devices on their own while enjoying enhanced sound clarity and extended battery life. The PerfectRite® works with a specific cleaning liquid cartridge & drying capsule, your patients will be able to buy refills at their preferred Audiologist Center.



Distributed by Microsonic, Inc.
Ambridge, PA, USA
Phone: 1 877 376 7139
www.safe-touch.com

IMPORTANT NOTE:

The unit requires 1 desiccant, 1 cleaning cartridge, and 1 adapter for operation.

770



SAFETOUCH™ PERFECTRITE® RIC CLEANING SYSTEM CONFIGURATIONS

- 770 SafeTouch™ PerfectRite® RIC Cleaning System
Contents: 1 Desiccant
- 770AS PerfectRite® RIC Cleaning System with Starter Kit
Contents: 1 desiccant, 1 Cleaning cartridge, 1 SMALL size adapter (for eardome sizes 4-6)
- 770AM PerfectRite® Cleaning System with Starter Kit
Contents: 1 desiccant, 1 Cleaning cartridge, 1 MEDIUM size adapter (for eardome sizes 7-8)
- 770AL PerfectRite® Cleaning System with Starter Kit
Contents: 1 desiccant, 1 Cleaning cartridge, 1 LARGE size adapter (for eardome sizes 9-12)

771



SAFETOUCH™ PERFECTRITE® ADAPTER FOR UNIVERSAL RIC/RITE EARDOMES

The PerfectRite® Adapters are designed to fit all size of receivers manufactured by major brands.

- 771S SafeTouch™ PerfectRite® RIC Adapter, Small (sizes 4-6)
- 771M SafeTouch™ PerfectRite® RIC Adapter, Medium (sizes 7-8)
- 771L SafeTouch™ PerfectRite® RIC Adapter, Large (sizes 9-12)

772



SAFETOUCH™ PERFECTRITE® CLEANING CARTRIDGES

The liquid cleaning solution is formulated to effectively dissolve earwax and disinfect the receiver. Removed earwax and dirt are collected in two designated chambers in cartridge. PerfectRite® Cleaning Cartridges are disposable. Approx. 25 cycles/cartridge

- 772 SafeTouch™ PerfectRite® Cleaning Refill, Box of 4 cartridges

773



SAFETOUCH™ PERFECTRITE® CUSTOM ADAPTER KIT

The custom adapter fabrication kit enables use of custom RIC/RITE earmolds as well as ITEs with PerfectRite®. Four custom adapters can be fabricated per kit.

- 773 SafeTouch™ PerfectRite® RIC Custom Adapter Fabrication Kit for 4

767S



SAFETOUCH™ PERFECTRITE® DESICCANTS


Special size orange color desiccant for PerfectRite®. Approx 25-30 cycles. The silica turns white at the end of life cycle. Individually seal-packed.

- 767S SafeTouch™ Desiccant Refill, Box of 4

Miscellaneous & Educational

MICROSONIC®

260



SWIM PLUG LANYARDS

Attaches to pre-drilled custom plugs to protect against loss or separation. Available in ten colors. (White, Black, Blue, Red, Orange, Yellow, Clear, Green, Pink, Purple)

#260 Swim Plug Lanyards

600



MINI BATTERY TESTER

Versatile battery tester/key holder slips into pocket or purse. Precision molded, it includes a compartment for two spare batteries. Tests size 675 and 13 size batteries.

#600 Mini Battery Tester

602



ACTIVAIR BUTTON CELL TESTER

A durable office battery tester that gives an immediate indication whether a battery is good or bad. Tests all types of button cells.

#602 Activair Button Cell Tester

603



DIGITAL HEARING AID BATTERY TESTER/HOLDER

Compact design is easy to carry, and gives an LCD readout of battery life remaining. Contains storage compartment for spares. For use with all sizes of hearing aid batteries.

#603 Digital Hearing Aid Battery Tester/Holder

608




ETY-PLUG HI-FI GENERIC ATTENUATORS

Designed to provide moderate attenuation in all sound levels. The ETY-PLUG has balanced attenuation which preserves the fidelity of the original sound. Noise reduction rating is 12db. One size fits all. Specify color: Beige, Translucent Red, Translucent Blue, Translucent Beige, Translucent Purple or Clear.

#608 ETY-PLUG Attenuators

#608A Cords for the ETY-PLUG Attenuators

609-610



EAR PLANES™

Disposable earplugs that reduce ear discomfort while flying. For travelers with colds, allergies or sinus conditions. Available in adult and children's sizes.

#609 Ear Planes™

#610 Children's Ear Planes™

611



COMMUNICATION COLLAR CLIP

Use with communication earmolds to clip tubing to collar and out of sight.

#611 Communication Collar Clip



TELEPHONE PADS™

Telephone pads reduce background noise and stop annoying feedback whistle during telephone conversation. Available in both round and square. Sold singly or in box of 24.


- #616 Telephone Pads™ (Round)
- #617 Telephone Pads™ (Round) (Box of 24)
- #618 Telephone Pads™ (Square)
- #619 Telephone Pads™ (Square) (Box of 24)



HEARING AID SWEAT BANDS

Protects BTE hearing aids from moisture and dust. Special fabric both absorbs and repels moisture but does not interfere with passage of sound. 5 Sizes(Mini,S,M,L,XL) 3 Mini Slim Sizes(Mini, S,M) (Specify size and color)

- #621 Hearing Aid Sweat Bands



SUPER SEALS™

A balloon-like protective sleeve for the BTE hearing aid that seals out most moisture and damaging elements such as body chemicals, hair spray, dyes, dirt and dust.

Available in 3 sizes – mini, standard and large
Available colors – neutral, beige, white, gray, natural, light brown, and brown
Please specify size and color when ordering.

- #625 Super Seal™ Kit (12 Seals & 1 Tool)
- #626 Refill package (12 Seals)
- #627 Expander Tool



EAR BAND-IT™

Swimmer's reversible neoprene headband with velcro-like closure in back to adjust for size. To be effective, it must be worn with earplugs or swim molds. Three sizes and five colors are available. Royal Blue, Red, Aqua, Pink, Purple.

- #631 Ear Band-It™ Small 1yr to 3 yrs
- #632 Ear Band-It™ Medium 4 yrs to 10 yrs
- #633 Ear Band-It™ Large 10 yrs to Adult



PERFECT EARMOLD IMPRESSIONS DVD

Learn about traditional and CIC impression techniques. In "Perfect Earmold Impressions" you'll see the entire impression making process, plus simple techniques for reducing remakes and increasing patient acceptance of earmolds. For all hearing health care professionals and educators. (22 minutes)

- #638 Perfect Earmold Impression DVD



MOUNTED DEMONSTRATION EAR

This beautiful silicone ear model is mounted on a handsome white pedestal. Suitable for display on desks, earmold tables and in offices. Available in two heights.

- #641 3-3/4" (95mm high) Mounted Demonstration Ear
- #642 6-1/2" (165mm high) Mounted Demonstration Ear



RESEALABLE PLASTIC BAGS (Package of 100)

Multi-purpose, resealable 3"x 5-1/2" plastic bags are handy for storage of small items. White labeling area on front.

- #643 Resealable Plastic Bags



INSTRUCTION EARS

These transparent model ears are ideal for instructing students in impression techniques. The canal is entirely visible so that proper placement of the oto-block or cotton block can be easily seen. One Pair: Available in white, beige and clear.

- #646 Instruction Ears



PREMIUM MAILING BOXES

3-1/2"x 2"x 2" hard cardboard box with wire clasps. Perfect for shipping earmold impressions. Sold individually. (Specify quantity)

- #647 Premium Mailing Boxes



SEE-THRU PLASTIC POUCHES

Colorful vinyl pouch, ideal for storing swim or sound plugs when not in use. Specify color (red, blue, green, yellow).

- #648 See-Thru Plastic Pouches



LEATHERETTE POUCH

Convenient pouch for storage of hearing aids, hearing protection molds, swim plugs and small parts.

- #659 Leatherette Pouch



MICROSONIC MAILING BOX

Postage-paid mailing box for sending impressions to Microsonic.

- #661 Microsonic Mailing Box



662

EARMOLD ORDER BOOK

Book of 50 Earmold Order Forms
#662 Earmold Order Book



663

663A

E-A-R™ CLASSIC EARPLUGS

These popular, disposable earplugs offer significant hearing protection when our custom Sound Defenders cannot be used.

#663 E-A-R™ Classic Earplugs Pkg. of 2
#663A E-A-R™ Classic Earplugs Box of 200 pairs



663B


663C

E-A-R™ EXPRESS POD PLUGS

These plugs add the push-in convenience of premolded plugs to the acoustic benefits of the traditional slow-recovery foam. Insertion and removal is quick and easy. Available with or without cords in boxes of 100 pairs. (NRR=25dbz)

Uncorded, box of 100 pairs, assorted colors Corded, box of 100 pairs, assorted colors.

#663B E-A-R™ Express Pod Plugs - uncorded pkg. of 2
#663C E-A-R™ Express Pod Plugs - uncorded box of 100 pairs
#663D E-A-R™ Express Pod Plugs - corded pkg. of 2
#663E E-A-R™ Express Pod Plugs - corded box of 100 pairs




670

671

GATOR CLIPS

Lightweight and nearly invisible, Gator Clips help reduce BTE hearing aid loss due to undetected slipping from the ear. Attach loop at one end to hearing aid, and the alligator clip at the other end to clothing. Available in single or bilateral styles.

#670 Gator Clip (*single*)
#671 Gator Clip (*bilateral*)



672

673

GATOR CLIPS II

Lightweight and nearly invisible, these reduce ITE or CIC hearing aid loss due to undetected slipping from the ear. Attach clasp at one end to hearing aid, and the alligator clip at the other end to clothing. Available in single or bilateral styles. Please specify.

#672 Gator Clips II (*single*)
#673 Gator Clips II (*bilateral*)




687

QUICK REFERENCE GUIDE TO EARMOLD SELECTION

Our convenient wall chart and accompanying technical bulletins will help you quickly select appropriate earmold styles.

#687 Quick Reference Guide to Earmold Selection



809

REPLACEMENT HANDLES

For easy removal of swim or sound plugs. Handles fit into holes pre-drilled at time of manufacture. Available in Red, Blue and Clear.

#809 Replacement Handles



817

UNIVERSAL HEARING SCREENER™

The screener is ideal for use in doctor's office and clinics. A LED lights when each frequency button is activated (500hz to 4000hz). Uses 2 AAA Batteries.

#817 Acrylic Hearing Aid Stands



900-902

901H

ER ATTENUATORS

Replacement attenuators for Musicians Earplugs. Available singly or in pairs, and in various colors as shown in chart. Available in beige, clear, brown, red, and blue.

#900 ER-25 Attenuators
#901 ER-15 Attenuators
#901H Solid Filters
#902 ER-9 Attenuators

	ER-25	ER-15	ER-15SP	ER-9
Beige	X	X		X
Brown	X	X		X
Clear	X	X		X
Orange			X	
Red	X	X		X
Blue	X	X		X



8009

8009A

SELF-INKING STAMP

Self-inking stamp with your account number, name and address.

#8009 Self-Inking Stamp
#8009A Black Ink



637

MICROSONIC CUSTOM EARMOLD MANUAL – 9TH EDITION

Revised in 2012 to cover most up-to-date information on earmold styles and materials. Microsonic Earmold Manual is widely recognized as the standard reference on the subject. It is used extensively in university programs as well as in dispensing offices. ISBN 978-1-4675-2187-1

#637 Microsonic Custom Earmold Manual – 9th Edition



684

685

STUDY GUIDE TO ACCOMPANY MICROSONIC'S CUSTOM EARMOLD MANUAL

Prepared specifically for the audiology student or other newcomer to hearing aid dispensing, this guide highlights the most important items in Microsonic's Custom Earmold Manual relating to good dispensing practice. It is a great review for the experienced practitioner, as well.

#684 Custom Earmold Manual Study Guide
#685 Custom Earmold Manual Study Guide - Instructor's Edition.

Microsonic
2960 Duss Avenue
Ambridge, PA 15003

Tel: 1-800-523-7672
Tel: 724-266-9270
Fax: 724-266-6309
www.microsonic-inc.com

Accessory Division
Tel: 1-877-376-7139
Fax: 724-266-9482
Online Store: <http://store.microsonic-inc.com>

