HEARING AID **PERFORMANCE / FEATURE** OVERVIEW

		mind440/PA	0440					
	40	10/04	4-775	.0	PA-11	0	4-705	3
	ARA	704×	9/08	7433	MA	⁷ 425	9/4/2	76
PLATFORM	CLE	mi	Int	Mi	4	Mi	Flag	RE
C-ISP	•							
Dual ISP		•		•		•		
ISP			•		•		•	•
CONNECTIVITY								
Wireless connectivity WidexLink	•							
InterEar (IE)	•							
CHANNELS								
	15	15	15	10	10	5	5	3
AUDIBILITY/ LISTENING IN QUIET								
IE TruSound Compression	•							
TruSound Compression		•		•				
EDRC	•	•	•	•	•	•	•	•
TruSound AOC	•	•		•				
TruSound Stabilisers	•	•		•				
Sound Stabiliser			•		•	•	•	•
Hammer Detector	•	•	•	•	•	•	•	•
TruSound Softener	•							
LISTENING IN NOISE								
High Definition Locator (multidirectional)	•	•	•	•	•	•		
Locator (Broadband)							•	•
Digital Pinna	•							
FreeFocus (via M-DEX)	•							
IE Speech Enhancer SII	•							
Speech Enhancer SII		•	•					
Noise reduction + SIS	•	•	•	•	•	•	•	
Noise reduction								•
Wind noise control	•							
ENHANCED USER COMFORT								
IE Feedback Cancelling	•							
Multidirectional Active Feedback Cancelling		•	•	•	•	•		
Feedback Cancelling							•	•
AISA	•	•	•	•	•	•		

	CLEAR440 mind440/PA-440 Inteo/PA-115 mind330 AIKIA/PA-110 mind220					?	\$2		
EXTENDED	CLEAR440	mind440/1	Inteo/PA-1	mind330	AIKIA/PA-	mind220	Flash/PA-10	REAL	
AUDIBILITY			_	-					
Audibility Extender	•	•	•	•		•			
ClearBand	•	•		•					
100 Hz - 10.5 kHz*		•		•					
70+ Hz - 11.2 kHz*	•								
USER CONTROLS				_					
IE Partner Alarm	•								
IE volume control	•								
IE program shift	•								
Listening programs	5	5	5	4	2+T/ MT	3	2+T/ MT	2+ T/MT	
Phone+	•								
IE Zen	•								
Zen		•		•		•			
Zen timer	•								
SmartSpeak	•	•		•		•			
SmartTone	•	•		•		•			
Acoustic beeps			•		•		•	•	
FITTING AND FINE TUNING									
Basic Sensogram	•	•	•	•	•	•	•	•	
Expanded Sensogram	•	•	•	•	250Hz	•			
In-situ vent effect	•	•	•	•	•	•			
In-situ RECD	•	•	•(m)	•(m)	•(m)				
QuickFit	•	•		•		•			
PrecisionFit	•	•	•	•	•	•			
SoundTracker	•	•	•	•	•	•	•	•	
Occlusion manager	•	•	•	•	•	•	•		
Sound Diary/Datalogging	•	•	•	•	•	•	•		
Event log	•	•	•	•	•	•			
PROGRAMMING									
NOAHlink		•	•	•	•	•	•	•	
nEARcom with NOAHlink	•								
USB Link	•								
iP5		•	•	•	•	•	•	•	
HI-PRO		•	•	•	•	•	•	•	

^{*} Frequency range is model dependent

REMOTE CONTROLS

RC3-1 BASIC

FEATURED

RC3-2 FULLY

RC4-1 BASIC

FEATURED

RC4-2 FULLY

COMPATIBLE WITH:	Inteo, AIKIA, PA-115 and PA-110	Inteo, AIKIA, PA-115 and PA-110	PA-440, mind440, mind330 and mind220	PA-440, mind440, mind330 and mind220
PROGRAM SELECTION	Choose from up to 3	Choose from up to 5	Choose from up to 3	Choose from up to 5
ON/OFF	-	•	•	-
FOR MONAURAL OR BINAURAL USE	•	•		•
VOLUME ADJUSTMENT	•	•	•	•
MUTE		•		•
ADJUSTABLE TIME-OUT	•	•	•	-
ALARM FUNCTION	•	•		•
EVENT-LOG FUNCTION		•	•	-
PROGRAM KEY	•			
SIMULTANEOUS ADJUSTMENT WITH BINAURAL FITTINGS	•		•	
SIMULTANEOUS OR SEPARATE ADJUSTMENT WITH BINAURAL FITTINGS		•		•
ZEN+ SMART TOGGLE*				-

^{*} In selected hearing aid models only

DEX - ASSISTIVE LISTENING DEVICES







TV-DEX



M-DEX

COMPATIBLE WITH CLEAR440

- Hearing aid remote control with LED (range 1 m)
- Volume up/down
- Program toggle
- Approx. 12 months of battery lifetime
- Key ring attachment

COMPATIBLE WITH CLEAR440

- Wireless device purpose-built for listen-ing to TV and audio. Range up to 10 m between Base and Controller.
- Up to 10 hrs of uninterrupted TV and stereo sound EchoFree audio quality
- Rechargeable in the TV-Base
- Volume up/down'Room Off'
- Transmission on/off

COMPATIBLE WITH CLEAR440

- · Wireless device for mobile phone applications, as well as advanced remote control with colour display
- On/off
- Accept/end call
- Navigation keys
 'Room Off'
- FreeFocus
- Integrated telecoil

