

Hearing Aid Direct Input Shoes & Leads

Most hearing aids have a direct connection facility, usually referred to as direct input. This facility allows the direct connection of the hearing aid to other audio equipment via a connector usually referred to as a shoe; thus the name 'Hearing aid direct input shoe'.


A direct connection into your hearing aid will give you the best possible sound quality and works brilliantly for listening to music. If you're new to direct input, we recommend you read our 'Beginners Guide to Direct Input' on the following page.

Direct input shoes are specific to the make and model of each hearing aid so it's important to get the right one.

We have grouped the shoes alphabetically by manufacturer. To find the correct shoe for your hearing aid; start by browsing through this section until you find the shoes relating to your hearing aid manufacturer and then look through the blue headings to find the correct model – keep in mind that many hearing aid models have very similar names, but require different shoes. For example, there are different shoes for the 'Phonak Naída UP' and the 'Phonak Naída SP'. If in doubt, please contact us for assistance before ordering.

We show all the shoes here in life size (as close as possible!) so you can compare an existing shoe against the page to help with identification.

Please note that although we have usually shown two shoes in the picture to give different views, the shoes are all sold as singles.



Battery size

By each shoe we indicate the size battery used by the compatible hearing aids. This is to help identification between different models of hearing aid with similar-sounding names.



Compatibility with direct input leads

By each shoe we indicate if it is suitable for use with direct input leads.



Suitable for use with direct input leads.



Not suitable for use with direct input leads.



Compatibility with ear level receivers

We also indicate whether it is suitable for use with 'ear level receivers' such as Roger X, MLxi, MicroMLxS & Amigo ear level receivers.



Suitable for use with ear level receivers.





Not suitable for use with ear level receivers.



Direct audio input leads

You will find direct input leads for use with your shoes at the end of this section – pages 64-65.

There are several different types of lead and the one that you require will depend on both your hearing instrument and what you want to listen to.

To listen to:	 fmGenie radio aid system	 iPods, tablets etc.
Hearing Instrument:	Any BTE hearing aid	All BTE hearing aids other than those made by Oticon
Solution:	fmGenie direct input leads	Blue plug direct input leads

Any Oticon BTE hearing aid	Cochlear implant or BAHA
fmGenie direct input leads	fmGenie filter leads or audio adaptor



Visit the website!

New direct input shoes are updated regularly plus some older shoes are still available for online ordering. We even have video clips showing how to fit the shoes and what to do with some of the extra bits and bobs and tools!



Go to www.connevans.info/sound-video to watch our video clips online.

Hearing aids not direct input compatible?



Some of the smallest hearing aids don't have a direct input facility but there are usually still other ways to connect in for listening to music.

You might have a loop program which can be used with an inductive neckloop or there might be a wireless streaming device available. Talk to your audiologist to find out the options available for your aids.

A Beginners Guide to Hearing Aid Direct Input

Many people who wear hearing aids every day are not aware that they could plug their hearing aids into a headset socket. The sound quality is fantastic and it is particularly good for listening to music. The method of connecting into a headphone socket is referred to as 'direct input' and uses a 'direct input' lead and hearing aid connection known as a 'hearing aid direct input shoe'.

The majority of behind-the-ear hearing aids, including almost all those in the current NHS range, have a direct input facility. We would suggest asking your audiologist if your hearing aid has a direct input facility and if it has been activated. If your direct input is not turned on, it will take your audiologist a matter of moments to activate the relevant program.

So how does hearing aid direct input work? As you will have noticed, your hearing aid doesn't have a socket on it to allow you to connect in a cable. This is because the hearing aid would need to be larger to fit a socket plus 'holes' in the hearing aid itself would attract dirt, dust and condensation, none of which are good for hearing aids. Instead, you use a 'shoe' as a connecting interface between the hearing aid and the cable. The shoe clips onto the bottom of your hearing aid and has a small three-pin socket on the bottom which allows you to plug in a direct input lead. The direct input lead can then be connected to your audio source (fmGenie, iPod, MP3, tablet, CD player etc).

Shoes: Direct input shoes are specific to the make and model of each hearing aid so it's important to get the right one!

In this section, we have grouped the shoes alphabetically by manufacturer. To find the correct shoe for your hearing aid; start by browsing through this section until you find the shoes relating to your hearing aid manufacturer and then look through the blue headings to find the correct model – keep in mind that many hearing aid models have very similar names, but require different shoes. For example, there are different shoes for the 'Phonak Naída UP' and the 'Phonak Naída SP'. If in doubt, please contact us for assistance before ordering.



You will see from this section that hearing aid direct input shoes come in all shapes and sizes, but what they all have in common is that they have a standard three-pin socket on the bottom.



Direct input leads: There are several different types of direct input lead, depending upon where and how they are to be used. The lead options are fully explained on page 64, however all these leads have the same three-pin plug at one end. Cleverly each of these three pins are a different size, meaning the lead will only plug into the shoe one way round, so you can't connect it incorrectly! We have leads available for people who have either one or two hearing aids, in a variety of different lengths and in black or white.

At the other end of the direct input lead is a standard 3.5mm plug – this will fit almost all headphone sockets. The leads are most popular for use with tablets and iPods, but they can also be used with devices such as computers, CD players, digital voice recorders, DAB radios, metal detectors...the list is endless!

If you want to connect to a product with a larger socket, you probably just need a simple adaptor (try our handy headphone extension/conversion kit MX100321 see page 215)



Direct input shoes and personal stereo direct input leads allow you to plug in anywhere that a pair of headphones can go.

If you have a Baha or a cochlear implant you can still take advantage of the benefits of direct input, you just need a different lead.

See pages 64-65 for more information.

Should you need any more information on direct input or audio shoes, the Connevens Customer Service Team will be pleased to help.



iPod & earphones



iPod & DPV lead

GN RESOUND Azure (AZ60), Danalogic6 (6060 & 4060), Metrix (MX60), Pixel (PL60), Plus5 (RP60)
D/I leads & ear level receivers


ID Marks: None
Contacts: Four
Colour: Translucent
Manfr. no: 15832600

DGNRA


Each	inc VAT	ex VAT
1+	£ 12.78	£ 10.65

GN RESOUND Danalogic6 (6090), Sparx
D/I leads & ear level receivers

ID Marks: None
Contacts: Four and a small battery contact pin in base, requiring the use of matching battery door (supplied)
Colour: Translucent
Manfr. no: 15349003

Please note that plugging in direct input leads requires a firm push (wiggling the plug may help).



Supplied with a replacement battery door unit (plus contact pins) which locks the shoe securely in place.


DGNRB


Each	inc VAT	ex VAT
1+	£ 12.78	£ 10.65

GN RESOUND Canta7 (780D), Canta2 (280D), Danalogic 283D
D/I leads & ear level receivers


ID Marks: GN ReSound DAI 11
Contacts: Four
Colour: Translucent
Manfr. no: 15072200

Please note that the hearing aid may require modification to the battery tray. If the shoe does not fit onto the hearing aid please contact your clinic to arrange modification.

DGNR780D


Each	inc VAT	ex VAT
1+	£ 12.78	£ 10.65

Technology moves fast – even more shoes available online

With the constant development of new hearing aids, direct input shoes are an ever-changing range. Our website is regularly updated with details of newly available shoes.

In addition to the models shown here, Connevens also has stocks of many shoes for non current hearing aids. If you can't find the shoe for the hearing aid you're looking for in the catalogue, try looking online.

Use the search box to find the correct shoe for your aid – enter one word of the make or model of your hearing aid.



www.DeafEquipment.co.uk

OTICON Sumo DM, Sumo XP, Sumo E, Spirit 3 Superpower, Swift 120+
D/I leads AP700


ID Marks: Oticon DAI AP-700
Contacts: Three
Colour: Grey
Manfr. no: 399-50-900-08

There is a small volume control which should be set to no. 4 when used with fmGenie equipment. Also has a HI-LO control which should be set to LO. This shoe attaches to the hearing aid by sliding onto the battery door.

DOAP700


Each	inc VAT	ex VAT
1+	£ 21.06	£ 17.55

Ear level receivers FM7


ID Marks: Oticon FM7
Contacts: Three
Colour: Grey
Manfr. no: 399-50-920-02

This shoe attaches to the hearing aid by sliding onto the battery door.

DOFM7


Each	inc VAT	ex VAT
1+	£ 21.06	£ 17.55

2 Hearing Aid Direct Input

OTICON GoPro, Safran, Spirit 3, Spirit 3 Direct, Spirit 3 VC & Spirit 3 Power, Syncro, Tego Pro, Tego

D/I leads AP800

ID Marks: Oticon AP800
Contacts: A fine contact array
Colour: Translucent with grey base
Manfr. no: 399-51-480-00

There is a small volume control which should be set to no. 4 when used with fmGenie equipment. Also has a HI-LO control which should be set to LO.

Does NOT fit Spirit 3 Superpower which requires DOAP700 shoe.

Supplied with adaptor tool



DOAP800

Each	inc VAT	ex VAT
1+	£ 26.36	£ 21.97



Ear level receivers FM8

ID Marks: Oticon FM8
Contacts: A fine contact array
Colour: Translucent with grey base
Manfr. no: 399-51-470-00

Does NOT fit Spirit 3 Superpower which requires DOFM7 shoe.

Supplied with adaptor tool



DOFM8

Each	inc VAT	ex VAT
1+	£ 21.84	£ 18.20



OTICON Agil, Acto, Chili, Epoq, Get, Hit, Ino, Safari, Sensei, Sensei Pro, Spirit Zest mini, Spirit Zest, Spirit Zest Power, Spirit Zest Super Power, Spirit Zest Plus, Spirit Synergy (not Mini), Vigo, Vigo Connect

D/I leads AP900

ID Marks: Oticon AP-900
Contacts: Five
Colour: Grey
Manfr. no: 399-50-520-08

This shoe attaches to the hearing aid by sliding onto the battery door.

Note – there is no Phonak MLx compatible shoe available for the Spirit Zest Mini.



DOAP900

Each	inc VAT	ex VAT
1+	£ 27.30	£ 22.75



Ear level receivers FM9

ID Marks: Oticon FM9
Contacts: Five
Colour: Translucent
Manfr. no: 399-50-590-00

This shoe attaches to the hearing aid by sliding onto the battery door.

Note – there is no Phonak ear level receiver compatible shoe available for the Spirit Zest Mini.

DOFM9

Each	inc VAT	ex VAT
1+	£ 26.21	£ 21.84



Micro BTE Tubes & Domes available!

A range of domes and micro tubes for BTE hearing aids without an earmould.

We are able to offer the Phonak & Oticon ranges of micro BTE tubes and domes.

These consumable items are supplied as manufacturer 'like for like' replacements, making ordering easy and speedy for you.



Go to www.DeafEquipment.co.uk and select **Micro BTE Tubes and domes** in the **Hearing Aid First Aid section**.

Or see page 130-131

PHONAK Milo UP, Milo Plus UP, Naída UP, Naída S UP, Naída Q UP, Nathos UP, Nathos S+ UP, Ok! Plus UP, Sky Q UP

D/I leads & ear level receivers AS10

ID Marks: None
Contacts: Three



Supplied with a replacement battery door unit which locks the shoe securely in place and fitting tool

D/I leads	Colour	ID Marks
DAS10BE	Beige	054-0022-01
DAS10BR	Brown	054-0022-03
DAS10BK	Black	054-0022-06
DAS10TA	Taupe	054-0022-35
DAS10CL	Clear	054-0022-13



Each	inc VAT	ex VAT
1+	£ 16.97	£ 14.14

PHONAK Milo SP, Milo Plus SP, Naída SP, Naída S SP, Naída Q SP, Nathos M & SP, Nathos S+ SP, Ok! Plus M & SP, Sky Q SP

D/I leads & ear level receivers AS11

ID Marks: None
Contacts: Three



More colours available online

Supplied with a replacement battery door unit which locks the shoe securely in place and fitting tool



DAS11BE	Beige	054-0115-01
DAS11BR	Brown	054-0115-03
DAS11BK	Black	054-0115-06
DAS11TA	Taupe	054-0115-35
DAS11CL	Clear	054-0115-13



Each	inc VAT	ex VAT
1+	£ 16.97	£ 14.14

PHONAK Ambra microP & SP, Cassia microP & SP, Dalia microP & SP, Nathos S+ M, Solana microP & SP

D/I leads & ear level receivers AS13

ID Marks: AS13
Contacts: Three



Supplied with a replacement battery door unit which locks the shoe securely in place and fitting tool



DAS13BE	Beige	054-0536-01
DAS13PA	Palladium Silver	054-0536-40
DAS13TA	Mocha Taupe	054-0536-F7
DAS13CH	Chocolate	054-0536-E5



Each	inc VAT	ex VAT
1+	£ 16.97	£ 14.14

PHONAK Bolero Q P & SP

D/I leads & ear level receivers AS16

ID Marks: AS16
Contacts: Three



More colours available online

Supplied with a replacement battery door unit which locks the shoe securely in place and fitting tool



DAS16BE	Beige	054-0119-01
DAS16P8	Velvet Black	054-0119-P8
DAS16CL	Clear	054-0119-13



Each	inc VAT	ex VAT
1+	£ 16.97	£ 14.14

PHONAK Certéna micro, Exélia micro & Exélia ART micro, Milo micro, Milo Plus micro, Nathos micro, Nios micro, Versáta micro

D/I leads & ear level receivers AS12

ID Marks: AS12
Contacts: Three

The aid may have a small sticker, covering the contacts inside the battery drawer, which will need to be removed for the shoe to work!



More colours available online

Supplied with a replacement battery door unit which locks the shoe securely in place and fitting tool



DAS12BE	Beige	054-0025-01
DAS12BR	Brown	054-0025-03
DAS12PA	Palladium Silver	054-0025-06
DAS12TA	Taupe	054-0025-35
DAS12CL	Clear	054-0025-13



Each	inc VAT	ex VAT
1+	£ 16.97	£ 14.14

PHONAK Ambra M H2O, Bolero Q-M13, Cassia M H2O, Dalia M H2O, Naída S CRT, Naída Q RIC, Nathos S+ micro, Nios S H2O, Sky Q-M13 & RIC, Solana M H2O

D/I leads & ear level receivers AS15


ID Marks: AS15 A
(updated version)
Contacts: Three

The aid may have a small sticker, covering the contacts inside the battery drawer, which will need to be removed for the shoe to work!



More colours available online

Supplied with a replacement battery door unit which locks the shoe securely in place and fitting tool



DAS15ABE	Beige	054-0682-01
DAS15ABK	Black	054-0682-06
DAS15AF3	Stirling Silver	054-0682-F3
DAS15AE5	Chocolate	054-0682-E5
DAS15ACL	Clear	054-0682-13



Each	inc VAT	ex VAT
1+	£ 16.97	£ 14.14

Wide range of Phonak colours available!

Most of the Phonak shoes are now available in a very wide range of colours to match the colour of the aids.

We don't have room for them here, but to view the full range for each shoe, please visit our website:

www.DeafEquipment.co.uk



2 Hearing Aid Direct Input

SIEMENS Prisma 2VC+, Signia

D/I leads & ear level receivers



ID Marks: F3 or 'for FM' is marked inside in the bottom.

Contacts: Four
Colour: Translucent
Manfr. no: 5630350

DSSIG

Each	inc VAT	ex VAT
1+	£ 16.20	£ 13.50

SIEMENS Infiniti Pro

D/I leads & ear level receivers

ID Marks: G5 inside shoe
Contacts: A fine array of 4 contacts
Colour: Translucent
Manfr. no: 7770584

This shoe is identical to DSIEO (you will notice it is ID'd G5), except for the colour coded contact covers for the hearing aid which are model specific.



Supplied with contact covers for the hearing aid.



DSIEB

Each	inc VAT	ex VAT
1+	£ 16.20	£ 13.50

SIEMENS Artis S, P & 2P, Centra HP, Chroma VC & P, Cielo P & 2P, Impact M & DP, Impact Pro M & DP, Prisma 2 Pro DP, Prisma 2 Pro P, Reflex DP

D/I leads & ear level receivers



ID Marks: L5 inside shoe
Contacts: Four
Colour: Translucent
Manfr. no: 10054830

DSIEG

Each	inc VAT	ex VAT
1+	£ 16.20	£ 13.50

SIEMENS Infiniti Pro SP, Music Pro SP, Prisma 2DSP, Prisma 2DSP+, Triano SP

D/I leads & ear level receivers



ID marks: T-SP or H5 inside
Contacts: Four
Colour: Translucent
Manfr. no: 7213296

DASTRI

Each	inc VAT	ex VAT
1+	£ 16.20	£ 13.50

SIEMENS Acuris S & P, Artis 2S, Artis 2SVC, Cielo 2, Centra P, S & SVC, Motion M & MVC, Triano S & SL

D/I leads & ear level receivers



ID marks: F5 inside shoe
Contacts: Four
Colour: Translucent
Manfr. no: 7220424

DSIED

Each	inc VAT	ex VAT
1+	£ 16.20	£ 13.50

SIEMENS Artis SP & 2SP, Centra SP, Chroma SP, Cielo 2SP, Impact & Impact Pro DSP, Intuis SP DIR, Nitro 300SP & 700SP, Prisma 2 Pro DSP & SP, Reflex DSP

D/I leads & ear level receivers



ID marks: M5 inside shoe
Contacts: Four
Colour: Translucent
Manfr. no: 10054617

DSIEI

Each	inc VAT	ex VAT
1+	£ 16.20	£ 13.50

Hearing aid batteries

Connevens stock a comprehensive range of hearing aid batteries from Rayovac, Duracell Activair & PowerOne, specially designed to be environmentally friendly. See Section 7 – page 122, or visit www.DeafEquipment.co.uk




SIEMENS Reflex M

D/I leads & ear level receivers

ID marks: G5
Contacts: A fine array of 4 contacts
Colour: Translucent
Manfr. no: 10299212

Please see note on DSIEB.



Supplied with contact covers for the hearing aid.



DSIEO




Each	inc VAT	ex VAT
1+	£ 16.20	£ 13.50

WIDEX Inteo IN-19, Senso Diva SD-19, Vita SV-19

D/I leads & ear level receivers Widex 5

ID Marks: Widex 5
Contacts: Five
Colour: Translucent
Manfr. no: 4 012 1900 0NA (295 W5)

With a slide on/off connection for the Widex SCOLA fm receiver



DW5



Each	inc VAT	ex VAT
1+	£ 20.26	£ 16.88

WIDEX Bravo B1/B2, Inteo IN-9E, Senso Diva SD9, Senso P7/P8, Vita SV-9

D/I leads & ear level receivers Widex 7

ID Marks: Widex 7 (previously Widex 1)
Contacts: 4, the centre left pin is offset and angled
Colour: Translucent
Manfr. no: 4 012 2700 050 (295 W7)

With a slide on/off connection for the Widex SCOLA fm receiver



DW7



Each	inc VAT	ex VAT
1+	£ 20.26	£ 16.88

SIEMENS Octiv HP & HP+, Octiv M & M+

D/I leads & ear level receivers

ID marks: P5
Contacts: Four



B G

DSIEQB Beige 10673021
DSIEQG Dark Granite 10674609



Each	inc VAT	ex VAT
1+	£ 18.72	£ 15.60

WIDEX Bravissimo BV-38, Bravo B32, Logo L32, Senso Vita SV-38, P37/P38,

D/I leads & ear level receivers Widex 6

ID Marks: Widex 6 (previously Widex 4, 3II, 3 and before that Logo Senso FM)
Contacts: 3 on one side, 1 on the other
Colour: Translucent
Manfr. no: 4 012 2400 0NA (295 DW6)

Has small catch on one side to lock shoe in place. With a slide on/off connection for the Widex SCOLA fm receiver



DW6



Each	inc VAT	ex VAT
1+	£ 20.26	£ 16.88

WIDEX Bravo B11/B12, Logo L6-L12 Logo Senso C series

D/I leads & ear level receivers Widex 8

ID Marks: Widex 8
Contacts: 3 on one side, 1 on the other
Colour: Translucent
Manfr. no: 4 012 3100 0NA (295 W8)

Fits Senso C8+ – C19+.

Has small catch on one side to lock shoe in place. With a slide on/off connection for the Widex SCOLA fm receiver



DW8



Each	inc VAT	ex VAT
1+	£ 20.26	£ 16.88

2 Hearing Aid Direct Input

Direct audio input leads

There are several different types of direct input lead, depending upon where and how they are to be used. Please use the chart below to select the correct lead:

To listen to:	 fmGenie radio aid system	 iPods, tablets etc.		
Hearing Instrument:	Any BTE hearing aid	All BTE hearing aids other than those made by Oticon	Any Oticon BTE hearing aid	Cochlear implant or BAHA
Solution:	fmGenie direct input leads	Blue plug direct input leads	fmGenie direct input leads	fmGenie filter leads or audio adaptor



fmGenie 3.5mm to Europlug single lead (1 aid)



fmGenie – Single lead for someone with 1 hearing aid



Length	400 mm		600 mm		800 mm		1 m	
Colour	Black or White		Black or White		Black or White		Black or White	
Part no.	FMGS400 FMGS400W		FMGS600 FMGS600W		FMGS800 FMGS800W		FMGS1M FMGS1MW	
	inc VAT	ex VAT	inc VAT	ex VAT	inc VAT	ex VAT	inc VAT	ex VAT
Qty 1-9	£ 13.75	£ 11.46	£ 14.05	£ 11.71	£ 14.38	£ 11.98	£ 14.69	£ 12.24
Qty 10-19	£ 12.37	£ 10.31	£ 12.65	£ 10.54	£ 12.94	£ 10.78	£ 13.22	£ 11.02
Qty 20-49	£ 11.69	£ 9.74	£ 11.94	£ 9.95	£ 12.22	£ 10.18	£ 12.48	£ 10.40

fmGenie 3.5mm to Europlug binaural lead (2 aids)



fmGenie – Binaural 'V' lead for someone with 2 hearing aids



Length	400 mm		600 mm		800 mm		1 m		1.5 m	
Colour	Black or White		Black or White		Black or White		Black or White		Black	
Part no.	FMGV400 FMGV400W		FMGV600 FMGV600W		FMGV800 FMGV800W		FMGV1M FMGV1MW		FMGV1M5	
	inc VAT	ex VAT	inc VAT	ex VAT	inc VAT	ex VAT	inc VAT	ex VAT	inc VAT	ex VAT
<i>Qty</i> 1-9	£ 22.16	£ 18.47	£ 22.73	£ 18.94	£ 23.28	£ 19.40	£ 23.81	£ 19.84	£ 25.48	£ 21.23
<i>Qty</i> 10-19	£ 19.94	£ 16.62	£ 20.46	£ 17.05	£ 20.95	£ 17.46	£ 21.43	£ 17.86	£ 22.93	£ 19.11
<i>Qty</i> 20-49	£ 18.84	£ 15.70	£ 19.32	£ 16.10	£ 19.79	£ 16.49	£ 20.23	£ 16.86	£ 21.66	£ 18.05



Blue plugged personal stereo single leads (1 aid)

– allows a single hearing aid to 'listen' to both channels of a stereo sound source.

Blue plug – Single lead for someone with 1 hearing aid



Length	800 mm		1 m	
Colour	Black or White		Black or White	
Part no.	DPS800 DPS800W		DPS1M DPS1MW	
	inc VAT	ex VAT	inc VAT	ex VAT
Qty 1	£ 23.88	£ 19.90	£ 24.50	£ 20.42
Qty 2-9	£ 21.49	£ 17.91	£ 22.06	£ 18.38
Qty 10-49	£ 20.30	£ 16.92	£ 20.83	£ 17.36

Blue plugged personal stereo V leads (2 aids)

– full stereo signal.



Blue plug – Binaural 'V' lead for someone with 2 hearing aids



Length	600 mm		800 mm		1 m		1.5 m	
Colour	Black or White		Black or White		Black or White		Black or White	
Part no.	DPV600 DPV600W		DPV800 DPV800W		DPV1M DPV1MW		DPV1M5 DPV1M5W	
	inc VAT	ex VAT	inc VAT	ex VAT	inc VAT	ex VAT	inc VAT	ex VAT
Qty 1	£ 26.34	£ 21.95	£ 26.90	£ 22.42	£ 27.53	£ 22.94	£ 30.49	£ 25.41
Qty 2-9	£ 23.71	£ 19.76	£ 24.22	£ 20.18	£ 24.78	£ 20.65	£ 27.44	£ 22.87
Qty 10-49	£ 22.39	£ 18.66	£ 22.87	£ 19.06	£ 23.40	£ 19.50	£ 25.92	£ 21.60

Direct input leads with inline RF interference filters for use with Cochlear implants and BAHAs

To listen to an fmGenie or to listen to music on battery powered audio devices such as a CD player or iPod, use these leads with the iConnect earhook or plug them directly into a Cochlear BP100, BP110 or Baha 4.



Advanced Bionics iConnect earhook

For use with AB Auria & Harmony cochlear implants. This special earhook attaches to the Auria/Harmony so that users can plug the fmGenie filter leads directly into their processor.

Adult size

3PAICE + colour code

Each	inc VAT	ex VAT
1+	£ 163.80	£ 136.50



Paediatric size

3PAICEP + colour code

Each	inc VAT	ex VAT
1+	£ 163.80	£ 136.50



fmGenie filter – Single lead for 1 hearing aid

	800 mm		1 m	
Length	800 mm		1 m	
Colour	Black		Black	
Part no.	FMGFS800		FMGFS1M	
	inc VAT	ex VAT	inc VAT	ex VAT
Qty 1-9	£ 20.18	£ 16.82	£ 20.51	£ 17.09
Qty 10-19	£ 18.17	£ 15.14	£ 18.46	£ 15.38
Qty 20+	£ 17.16	£ 14.30	£ 17.44	£ 14.53

fmGenie filter – Binaural 'V' lead for someone with 2 hearing aids

	800 mm		1 m		1.5 m	
Length	800 mm		1 m		1.5 m	
Colour	Black		Black		Black	
Part no.	FMGFV800		FMGFV1M		FMGFV1M5	
	inc VAT	ex VAT	inc VAT	ex VAT	inc VAT	ex VAT
Qty 1-9	£ 33.35	£ 27.79	£ 33.90	£ 28.25	£ 35.54	£ 29.62
Qty 10-19	£ 30.01	£ 25.01	£ 30.52	£ 25.43	£ 31.99	£ 26.66
Qty 20+	£ 28.34	£ 23.62	£ 28.81	£ 24.01	£ 30.22	£ 25.18

Cochlear Implant personal audio cables

Allow you to connect a personal audio device directly into your CI processor. Please note these are only suitable for use with battery powered equipment such as iPods and personal CD players because they do not have an isolating transformer in them to protect you.



Cochlear personal audio cables

For use with	Cochlear Freedom		Cochlear Nucleus 5 & 6	
Type	Monaural		Monaural or Binaural	
Part no.	DCPACFR		DCPACCP800 DCPACCP800B	
	inc VAT	ex VAT	inc VAT	ex VAT
Qty 1+	£ 55.73	£ 46.44	£ 66.10	£ 55.08

BAHA audio adaptors

Allow you to connect a personal audio device directly into your sound processor. The audio adaptor has an inbuilt isolation transformer which protects the wearer when connecting the sound processor to external sound equipment.



BAHA audio adaptors

For use with	Cochlear Baha Divino, Intenso, Classic 300 & Compact		Cochlear Baha Cordell II		Cochlear Baha BP100, BP110 & Baha4		Oticon Medical Ponto, Ponto Pro, Ponto Pro Power	
Type	Monaural		Monaural		Monaural		Monaural	
Length	1.3 m		1.9 m		0.7 m		1.3 m	
Part no.	DCBA065		DCBA067		DCBA91180		DCOM50377	
	inc VAT	ex VAT	inc VAT	ex VAT	inc VAT	ex VAT	inc VAT	ex VAT
Qty 1+	£ 132.19	£ 110.16	£ 132.19	£ 110.16	£ 132.19	£ 110.16	£ 166.90	£ 139.08

Magnifying Glasses – ideal for teachers working with many different shoes

Fitting and identifying direct input shoes can be a fiddly job. Markings are often small and hard to read. We have an excellent range of magnifying glasses available in different shapes and sizes!

Please visit www.DeafEquipment.co.uk and search 'magnify' or see page 135 for more details.

